

OPA

High-level disinfectant for medical devices, heat-sensitive dental & endoscopic instruments


HISEPT OPA is a high-level disinfectant for reusable medical devices, patient-care equipment, heat-sensitive dental and endoscopic instruments. This product contains ortho-phthalaldehyde 0.575% as a disinfection agent, with fast-acting and broad-spectrum antimicrobial activity. It eliminates virtually all recognized pathogenic microorganisms on medical devices, including all steam and heat-sensitive semi-critical instruments.

Characteristics

- Contains 0.575% Ortho-Phthalaldehyde
- Broad-spectrum antimicrobial activity
- Ready to use without the need of any activator
- Reusable for up to 14 days

Product Specification

Color/Form	Clear & colorless solution
pH	6 - 8
Scent	Characteristic odor of aldehyde
Shelf Life	24 months



THIS IS A HISEPT PRODUCT, WITH
**HIGH-LEVEL
SAFETY &
PROTECTION**

Ingredients

- 100 g HISEPT OPA contains:
- 0.575 g ortho-phthalaldehyde

Auxiliary ingredients:

- Chelating agents
- Buffering agent
- Corrosion inhibitors
- Deionized water

Microbiological Effectiveness

Meets European Standards

Bactericidal Activity

According to EN 13727 (contact time: 10 minutes)	Staphylococcus. aureus/ ATCC 6538 Enterococcus. hirae/ ATCC 10541 Pseudomonas. aeruginosa/ ATCC 15442 Escherichia. coli/ ATCC 10536
According to EN 14561 (contact time: 10 minutes)	Staphylococcus. aureus/ ATCC 6538 Enterococcus. hirae/ ATCC 10541 Pseudomonas. aeruginosa/ ATCC 15442

Yeasticidal Activity

According to EN 13624 (contact time: 10 minutes)	Candida. albicans/ ATCC 10231
According to EN 14562 (contact time: 10 minutes)	Candida. albicans/ ATCC 10231

Fungicidal Activity

According to EN 13624 (contact time: 10 minutes)	Aspergillus. niger/ ATCC 16404
According to EN 14562 (contact time: 15 minutes)	Aspergillus. niger/ ATCC 16404

Sporicidal Activity

According to EN 13704 (contact time: 60 minutes)	Bacillus. subtilis/ ATCC 6633
--	-------------------------------



Directions for use

Immerse clean device completely into the ready-to-use HISEPT OPA, filling all lumens and eliminating bubbles, then cover the tank for at least 10 minutes. Remove device from the solution and rinse thoroughly with clean water.

Compatibility


HISEPT OPA is compatible with most materials such as stainless steel, copper, polyethylene, rubber, PVC, polyurethane and silicon.

Packagings

3 canisters per carton,	5 L each
6 bottles per carton,	1 L each

Hazard Rating

Ready-to-use	NFPA
Health	2
Flammability	0
Reactivity	0



Warning

- The disinfection tank shall always be covered with a lid. If the solution is kept in a closed container, it will be reusable for up to 14 days.
- In case of contact with skin and eyes, rinse immediately with plenty of water and seek medical advice.
- Keep out of reach of children and unaware persons.
- Store between + 5 °C and + 35 °C and away from direct sunlight
- Contact with Ortho-Phthalaldehyde solution may stain exposed skin or clothing. It is recommended to use nitrile or butyl rubber gloves. When using latex gloves, the user should use double gloves or change single gloves frequently.



Certifications

- >> ISO 9001: 2015
- >> ISO 13485: 2016
- >> GMP
- >> FSC



Middle East Office
Behban Pharmed Lotus Co.
Address : Unit 113, 12 Floor, Bldg. B, Kaveh
Complex, Azadi Ave., Tehran, Iran
Tel: (+9821) 66581294-5/66581269

info@hisept.com

World Locations: Düsseldorf, Kuala Lumpur, Istanbul, Tehran, Baghdad