

# SEPTOSTERILE-A

## High-level disinfectant for automated disinfection of flexible endoscopes and heat-sensitive medical devices


HISEPT SEPTOSTERILE-A is a high-level disinfectant for chemical reprocessing of flexible endoscopes & medical devices. This product contains glutaraldehyde 2.4% as a disinfection agent, with fast-acting and broad-spectrum antimicrobial activity. It eliminates virtually all recognized pathogenic microorganisms on medical devices, including all steam and heat-sensitive semi-critical instruments. HISEPT SEPTOSTERILE-A has good material compatibility especially with plastics and metallic materials.

## Characteristics

- Good material compatibility
- Broad-spectrum antimicrobial activity
- Ready-to-use without the need of any activator
- Limited evaporation of glutaraldehyde

## Product Specification

Color/Form	Clear green solution
pH	5.5 - 6.6
Scent	Characteristic odor of aldehyde and lemon scent
Shelf Life	24 months



THIS IS A HISEPT PRODUCT, WITH  
**HIGH-LEVEL  
SAFETY &  
PROTECTION**

## Ingredients

100 g HISEPT SEPTOSTERILE-A contains:

- 2.4 g glutaraldehyde

Auxiliary ingredients:

- Low-foam surfactant
- Buffering agent
- Corrosion inhibitors
- Essence
- Dye
- Deionized water

## Microbiological Effectiveness

### Meets European Standards

#### Bactericidal Activity

According to EN 13727 Staphylococcus. aureus/ ATCC 6538  
(contact time: 10 minutes) Enterococcus. hirae/ ATCC 10541  
Pseudomonas. aeruginosa/ ATCC 15442  
Escherichia. coli/ ATCC 10536

According to EN 14561 Staphylococcus. aureus/ ATCC 6538  
(contact time: 10 minutes) Enterococcus. hirae/ ATCC 10541  
Pseudomonas. aeruginosa/ ATCC 15442

#### Yeasticidal Activity

According to EN 13624 Candida. albicans/ ATCC 10231  
(contact time: 10 minutes)

According to EN 14562 Candida. albicans/ ATCC 10231  
(contact time: 10 minutes)

#### Fungicidal Activity

According to EN 13624 Aspergillus. niger/ ATCC 16404  
(contact time: 10 minutes)

According to EN 14562 Aspergillus. niger/ ATCC 16404  
(contact time: 30 minutes)

#### Sporicidal Activity

According to EN 13704 Bacillus. subtilis/ ATCC 6633  
(contact time: 60 minutes)

#### Virucidal Activity

According to EN 14476 Herpes simplex virus(HSV-1)  
(contact time: 10 minutes)

## Directions for use

**Cleaning:** The cleaning step plays a key role in chemical processing. The success of the processing depends on thorough cleaning in order that the flexible endoscope is free of secretion and protein residuals. This prevents the incrustation of any remaining contamination during the subsequent disinfecting process.

**Reprocessing of flexible endoscopes:** Please refer to the manufacturer's recommendations.

## Compatibility

HISEPT SEPTOSTERILE-A is compatible with most materials such as stainless steel, copper, plastic, rubber and silicon.

2/2 Latest revision: Dec 1, 2023

## Packagings

3 canisters per carton,


5 Leach

6 bottles per carton,

1 L each

## Hazard Rating

Ready-to-use	NFPA
Health	2
Flammability	0
Reactivity	0



## Warning

- The disinfection tank shall always be covered with a lid.
- In case of contact with skin and eyes, rinse immediately with plenty of water and seek medical advice.
- Keep out of reach of children and unaware persons.
- Store between + 5 °C and + 35 °C and away from direct sunlight.



## Certifications

- >> ISO 9001: 2015
- >> ISO 13485: 2016
- >> GMP
- >> FSC



**HISEPT**<sup>®</sup>  
HIGH-LEVEL SAFETY & PROTECTION

Middle East Office  
**Behban Pharmed Lotus Co.**  
Address : Unit 113, 12 Floor, Bldg. B, Kaveh  
Complex, Azadi Ave., Tehran, Iran  
Tel: (+9821) 66581294-5/66581269

info@hisept.com

**World Locations:** Düsseldorf, Kuala Lumpur, Istanbul, Tehran, Baghdad